Jazzy 1121

UK

AGN# 088 609 661







Your new Pride Illustrated Parts Breakdown (IPB) Manual has been compiled from the latest product information available at the time of publication. Pride Mobility Products® Corporation reserves the right to make changes and updates to this manual and to the product(s) it represents. Any changes, updates, or improvements to the product(s) may result in slight differences between this manual and the product you purchased.

An IPB is an illustrated parts listing of all product components which have maintenance significance. IPB's utilize "exploded" drawings to detail the individual components for each assembly and an itemized list of reference numbers, part numbers, descriptions, and quantities per assembly. Some parts may be included in an assembly but may not be available for purchase separately. Refer to the parts list to locate the part numbers for assemblies and components for ordering what you need.

To achieve proper identification and relationship of detail parts and next higher assemblies, parts listings are presented in a convenient alphanumeric reference system. Reference numbers for parent assemblies identify one part number that can be used to order all of the components listed within it. Components of an assembly will have the same alpha character as its parent assembly.

■ For example, the parent assembly part number using reference number of A1 — A4 is the only number you need to place an order for parts with references of A1, A2, A3, and A4. Instead of ordering 4 separate parts, use the parent assembly part number to get the complete assembly.

Master assembly part numbers have reference numbers that include several subassemblies. By ordering parts using master assembly part numbers, you can save a lot of time and reduce the chance for error.

■ For example, A1 — G10, a master assembly reference number, includes all components of assemblies A, B, C, D, E, F, and G. Use caution when ordering assemblies. Check to see if a master assembly part number exists for the components that you need. If an entire master assembly is not needed, then the individual assembly part number (i.e. for reference number G1 — G10) will be sufficient to use for ordering.

Asterisks (*) and NOTES are used throughout this service manual to point out important information regarding assembly structure and contents.

■ For example, in some illustrations, only the left OR right side is shown. In other drawings, only the front OR back is depicted. The part numbers for ALL components are listed. An asterisk will also be used to indicate which components are not illustrated. Illustrated Items without callouts are shown for assembly purposes.

Updates to the contents of this IPB will be made on a regular basis. Providers may retrieve updates immediately upon their release by going to the Pride Service Website. The web address is www.prideservice.com. We welcome any questions or comments regarding the new Pride IPB Manual that you may have.

Please note that all pricing that is included with this manual is subject to change without any prior notice from Pride.

Thank you for choosing products from Pride Mobility Products[®] Corporation.


Copyright © 2006 Pride Mobility Products® Corporation www.pridemobility.com INFMANU2112/Rev.O/120806

Table of Contents



SECTION I.	
Main Frame Assembly	4
SECTION II.	
Single Brake Release Assembly	6
SECTION III.	
Motor Assembly	8
Single Brake Release Lever, Standard	
Single Brake Release Lever, HD Option	
Individual Brake Release Lever, Standard	12
Individual Brake Release Lever, HD Option	14
SECTION IV.	
Anti-tip Assembly	
Standard	
Heavy Duty Option	
SECTION V.	
Articulating Beam Assembly	20
Straight	
Curved	
SECTION VI.	
Electronic Tray Assembly	24
VSI	
Remote Plus	
SECTION VII.	
Controller Assembly	28
SECTION VIII.	
Foot Platform Assembly	30
9.55"	
11.62	
Heavy Duty	
SECTION IX.	
Wheel Assembly	
Flat-free	
Pneumatic	37
SECTION X.	
Shroud Assembly	
Standard	
Power Elevating Seat	39
SECTION XI.	
Off-board Charger	40
SECTION XII.	
Owner's and Decal Packages	41
APPENDIX A.	
Electronic System Diagrams	42
VSI	
Remote Plus	43





I. Main Frame Assembly




REF#	PART#	DESCRIPTION	QTY
A1 — A15	FRMASMB4093	FRAME,ASSY,MAIN,JAZZY1101	1
A1	DWR1101E005	FRAME,FRAME WELDMENT,MAIN, JAZZY1101/1121,E-1101-005	1
A2	DWR1101C033	BUMPER,SKID BAR,MODIFIED FOR RELEASE HANDLE,C-1101-033	1
A3	DWR1101A025	BUSHING,STEEL,0.335" ID X 0.51" OD X 0.670"TH,A-1101-025	1
A4	DWR1101E012	FRAME,BATTERY DOOR WELDMENT,JAZZY1101/1121,E-1101-012	1
A5	DWR1101E019	BATTERY,DOOR,JAZZY1101/1121,E-1101-019	1
A6	DECSERN1032	DECAL,SERIAL NUMBER,4" X 6",BARCODE,W/PRIDE LOGO	1
*A7	PLSTIE1004	TIE WRAP,BLACK,8"	8
A8	DWR1101D018	BATTERY,SPACER,JAZZY1101/1121,D-1101-018	1
A9	CLISTL1012	CLIP,SNAP-IN,0.228"D	4
A10	DWR1101E022	SHROUD,BELLY PAN,RIGHT, JAZZY1101/1121,E-1101-022	1
A11	DWR1101E023	SHROUD,BELLY PAN,LEFT, JAZZY1101/1121,E-1101-023	1
A12	SCRPPH1050	SCREW,PPH,BLACK ZINC,M47 X 10MM	10
A13	SCRASBH1091	SCREW,BHSCS,BLACK ZINC,M8-1.25 X 140MM,GRADE 10.9	2
A14	SCRSHFH1041	SCREW,SHFH,BLACK ZINC,M8-1.25 X 55MM,GRADE 10.9	1
A15	SCRSETS1102	SCREW, SET,M6-1.0 X 6MM,CUP POINT,W/NYLON PATCH	4

Notes:

*A7 Not shown.








REF#	PART #	DESCRIPTION	QTY
A1 — A15	FRM ASM B4106	FREEWHEEL, ASSY, SINGLE BRAKE RELEASE LEVER, JAZZY1101/1121	1
A1	DWR1101D007	FREEWHEEL,HANDLE BRACKET WELDMENT,ANGLED, JAZZY1101/1121,D-1101-007	1
A2	DWR1101E017	FREEWHEEL, LEVER CABLE, 27.38"L, JAZZY 1101/1121, E-1101-017	2
A3	DWR1122B046	LATCH,CATCH LATCH,B-1122-046	1
A4	DWR1122A021	SPRING, DOUBLE TORSION, STAINLESS STEEL, 0.035", #302 WIRE, A-1122-021	1
A5	DWG1000A800	FREEWHEEL,SMALL L BRACKET,SINGLE LEVER BRAKE RELEASE STYLE #2,A-1000-080	2
A6	NUTNJAM1014	NUT,NYLOCK JAM,BLACK ZINC,1/4-20	6
A7	DWR1101A016	CABLE BARREL,BRIGHT ZINC,M47 TAP,A-1101-016	2
A8	BLTASHL1084	BOLT,SHOULDER,6MM X 10MM	2
A9	SCRASBH1071	SCREW,BHSCS,BLACK ZINC,M6-1.0 X 25MM,GRADE 8.8	1
A10	NUTNYLK1019	NUT,NYLOCK,BLACK ZINC,M6-1.0	1
A11	SCRPPH1050	SCREW,PPH,BLACK ZINC,M47 X 10MM	2
A12	DWG1200B770	FREEWHEEL,RELEASE LEVER,CABLE READY,B-1200-077	2
A13	DWR1101C008	FREEWHEEL,L BRACKET,JAZZY1101/1121,C-1101-008	1
A14	SCRSETS1026	SCREW,SET,M6-1.0 X 6MM	2
A15	PLSTIE1001	TIE WRAP,BLACK,4"	2

Available on units assembled prior to December 30, 2003.





III. Motor Assembly — Single Brake Release Lever, Standard




REF#	PART#	DESCRIPTION	QTY
A1 — B7	DRVASMB1513	MOTOR,ASSY,DRIVE,E620HS, MWD, W/SINGLE BRAKE, PAIR,MTS	1
A1a — B7	DRVASMB1438	MOTOR,ASSY,DRIVE,E620HS,MWD,RIGHT,MTS	1
A1b — B7	DRVASMB1437	MOTOR,ASSY,DRIVE,E620HS,MWD,LEFT,MTS	1
A1a	DWR1101D025	MOUNTING BRACKET,MOTOR MOUNT BRACKET WELDMENT,RIGHT, JAZZY1101/1121,D-1101-025	1
A1b	DWR1101D024	MOUNTING BRACKET,MOTOR MOUNT BRACKET WELDMENT,LEFT, JAZZY1101/1121,D-1101-024	1
A2	NUTNYLK1041	NUT,NYLOCK,ZINC,M10-1.25,GRADE 8	1
A3	WASSTND1138	WASHER,FLAT,BLACK ZINC,10.3MM ID X 24MM OD X 2.8MM TH,SERRATED,B-1005-033	1
A4	DWR1101A019B	BEARING, BALL,12.0MM ID X 28.0MM OD X 8.0MM TH,A-1101-019,(SKF BEARING: 6001-2RS DOUBLE SEAL)	2
B1a — B7	DRVMOTR1173	MOTOR,DRIVE,E620,EXTERNAL THREAD, NO HANDLE/SUSPENSION, RIGHT,H-1101-004,(MTS:62020023)	1
B1b — B7	DRVMOTR1172	MOTOR,DRIVE,E620,EXTERNAL THREAD, NO HANDLE/SUSPENSION, LEFT,H-1101-003,(MTS:62010023)	1
B1	INDPART1925	AXLE KEY,0.236"W X 0.236" H X 1.566" L	1
B2	PLSINDV1104	BRAKE,CAP,E610	1
В3	INDPART2056	BRAKE,SOLENOID,E610	1
B4	DRVBRSH1006	MOTOR,BRUSH,(MTS: 62400001)	2
B5	PLSINDV1103	MOTOR,BRUSH CAP,(MTS: 62400013)	2
В6	SCRALNH1037	SCREW,SHCS,M6-1.0 X 16MM	3
B7	SCRALNH1038	SCREW,SHCS,M6-1.0 X 30MM	3

Notes:

Available on units assembled prior to December 30, 2003. Top-level assemblies will include components designated with corresponding lower-case letters. Parts without lower-case extensions will be included in all top-level assemblies. This assembly is for use with single lever brake release assembly.





III. Motor Assembly — Single Brake Release Lever, HD Option




REF#	PART#	DESCRIPTION	QTY
A1 — B7	DRVASMB1537	MOTOR,ASSY,DRIVE,E610,W/ SINGLE BRAKE, MWD, PAIR,MTS	1
A1a — B7	DRVASMB1539	MOTOR,ASSY,DRIVE,E610,W/ SINGLE BRAKE, MWD, RIGHT,MTS	1
A1b — B7	DRVASMB1538	MOTOR,ASSY,DRIVE,E610,W/ SINGLE BRAKE, MWD, LEFT,MTS	1
A1a	DWR1101D025	MOUNTING BRACKET, MOTOR MOUNT BRACKET WELDMENT, RIGHT, JAZZY1101/1121,D-1101-025	1
A1b	DWR1101D024	MOUNTING BRACKET,MOTOR MOUNT BRACKET WELDMENT,LEFT, JAZZY1101/1121,D-1101-024	1
A2	NUTNYLK1041	NUT,NYLOCK,ZINC,M10-1.25,GRADE 8	1
A3	WASSTND1138	WASHER,FLAT,BLACK ZINC,10.3MM ID X 24MM OD X 2.8MM TH,SERRATED,B-1005-033	1
A4	DWR1101A019	BEARING, BALL, 12.0MM ID X 28.0MM OD X 8.0MM TH, A-1101-019, (SKF BEARING: 6001-ZZ DOUBLE SHIELD)	2
B1a — B7	DRVMOTR1205	MOTOR,DRIVE,E610,3700 RPM,RIGHT,MIP,H-1121-002, (MTS-62020078)	1
B1b — B7	DRVMOTR1204	MOTOR,DRIVE,E610,3700 RPM,LEFT,MIP,H-1121-003, (MTS-62010078)	1
B1	INDPART1925	AXLE KEY,0.236"W X 0.236" H X 1.566" L	1
B2	PLSINDV1104	BRAKE,CAP,E610	1
В3	INDPART2056	BRAKE,SOLENOID,E610	1
B4	DRVBRSH1006	MOTOR,BRUSH,(MTS: 62400001)	2
B5	PLSINDV1103	MOTOR,BRUSH CAP,(MTS: 62400013)	2
B6	SCRALNH1037	SCREW,SHCS,M6-1.0 X 16MM	3
B7	SCRALNH1038	SCREW,SHCS,M6-1.0 X 30MM	3

Notes:

Available on units assembled prior to December 30, 2003. Only for use with Jazzy1121HD package. Top-level assemblies will include components designated with corresponding lower-case letters. Parts without lower-case extensions will be included in all top-level assemblies.





III. Motor Assembly — Individual Brake Release Levers, Standard




REF#	PART#	DESCRIPTION	QTY
A1 — C11	DRVASMB1514	MOTOR,ASSY,DRIVE,E620HS, MWD, W/INDIVIDUAL BRAKE LEVERS, PAIR,MTS	1
A1a — C11	DRVASMB1516	MOTOR,ASSY,DRIVE,E620HS, MWD, INDIVIDUAL BRAKE LEVER, RIGHT,MTS	1
A1b — C11	DRVASMB1515	MOTOR,ASSY,DRIVE,E620HS, MWD, INDIVIDUAL BRAKE LEVER, LEFT,MTS	1
A1a	DWR1101D025	MOUNTING BRACKET,MOTOR MOUNT BRACKET WELDMENT,RIGHT, JAZZY1101/1121,D-1101-025	1
A1b	DWR1101D024	MOUNTING BRACKET,MOTOR MOUNT BRACKET WELDMENT,LEFT, JAZZY1101/1121,D-1101-024	1
A2	NUTNYLK1041	NUT,NYLOCK,ZINC,M10-1.25,GRADE 8	1
А3	WASSTND1138	WASHER,FLAT,BLACK ZINC,10.3MM ID X 24MM OD X 2.8MM TH,SERRATED,B-1005-033	1
A4	DWR1101A019	BEARING, BALL, 12.0MM ID X 28.0MM OD X 8.0MM TH, A-1101-019, (SKF BEARING: 6001-ZZ DOUBLE SHIELD)	2
B1a — B7	DRVMOTR1173	MOTOR,DRIVE,E620,EXTERNAL THREAD, NO HANDLE/SUSPENSION, RIGHT,H-1101-004,(MTS:62020023)	1
B1b — B7	DRVMOTR1172	MOTOR,DRIVE,E620,EXTERNAL THREAD, NO HANDLE/SUSPENSION, LEFT,H-1101-003,(MTS:62010023)	1
B1	INDPART1925	AXLE KEY,0.236"W X 0.236" H X 1.566" L	1
B2	PLSINDV1104	BRAKE,CAP,E610	1
В3	INDPART2056	BRAKE,SOLENOID,E610	1
B4	SCRALNH1037	SCREW,SHCS,M6-1.0 X 16MM	3
B5	SCRALNH1038	SCREW,SHCS,M6-1.0 X 30MM	3
B6	DRVBRSH1006	MOTOR,BRUSH,(MTS: 62400001)	2
B7	PLSINDV1103	MOTOR,BRUSH CAP,(MTS: 62400013)	2
C1a — C11	FRMASMB5479	FREEWHEEL,ASSY,RELEASE LEVER,JAZZY1101/1121,RIGHT	1
C1b — C11	FRMASMB5478	FREEWHEEL,ASSY,RELEASE LEVER,JAZZY1101/1121,LEFT	1
C1a	DWR1101D030	FREEWHEEL LEVER W/OVERMOLD COVER, RIGHT, D-1101-030	1
C1b	DWR1101D031	FREEWHEEL,LEVER W/OVERMOLD COVER,LEFT	1
C2a	DWR1133C014	FREEWHEEL,RELEASE ROD,RIGHT, JAZZY1133,C-1133-014	1
C2b	DWR1133C015	FREEWHEEL,RELEASE ROD,LEFT, JAZZY1133,C-1133-015	1
C3	DWR1114B014	FREEWHEEL,MOTOR RELEASE WELDMENT,JAZZY1113ATS,B-1114-014	1
C4	SCRSETS1023	SCREW,SET,M6-1.0T X 8MM,CUP POINT	1
C5	SPGRAW1046	SPRING,8MM ID X 10MM OD X 10MM TH,COMPRESSION,FREEWHEEL RELEASE LEVER,MTS	1
C6	WASWAVY1008	WASHER,BELLEVILLE,STAINLESS STEEL,0.317"ID X 0.625"OD X .022"TH	1
C7	PINSOTR1000	PIN, COTTER,ZINC PLATED,1/16" X 1"	2
C8	BLTASHL1111	BOLT,SHOULDER,8MM X 45MM,GRADE 12.9,C-1056-025	1
C9	NUTNYLK1019	NUT,NYLOCK,BLACK ZINC,M6-1.0	1
C10	WASNYLN1089	WASHER,NYLON,8MM ID X 19MM OD X 1MM TH	1
C11	WASSTND1100	WASHER,FLAT,ZINC,6MM ID X 14MM OD X 1.5MM TH	1

Notes

Top-level assemblies will include components designated with corresponding lower-case letters. Parts without lower-case extensions will be included in all top-level assemblies.





III. Motor Assembly — Individual Brake Release Levers, HD Option




REF#	PART#	DESCRIPTION	QTY
A1 — C11	DRVASMB1579	MOTOR,ASSY,DRIVE,E610,W/INDIVIDUAL BRAKE LEVERS, PAIR, MWD,MTS	1
A1a — C11	DRVASMB1581	MOTOR,ASSY,DRIVE,E610,W/INDIVIDUAL BRAKE LEVERS, RIGHT, MWD,MTS	1
A1b — C11	DRVASMB1580	MOTOR,ASSY,DRIVE,E610,W/INDIVIDUAL BRAKE LEVERS, LEFT, MWD,MTS	1
A1a	DWR1101D025	MOUNTING BRACKET,MOTOR MOUNT BRACKET WELDMENT,RIGHT, JAZZY1101/1121,D-1101-025	1
A1b	DWR1101D024	MOUNTING BRACKET,MOTOR MOUNT BRACKET WELDMENT,LEFT, JAZZY1101/1121,D-1101-024	1
A2	NUTNYLK1041	NUT,NYLOCK,ZINC,M10-1.25,GRADE 8	1
A3	WASSTND1138	WASHER,FLAT,BLACK ZINC,10.3MM ID X 24MM OD X 2.8MM TH,SERRATED,B-1005-033	1
A4	DWR1101A019	BEARING, BALL, 12.0MM ID X 28.0MM OD X 8.0MM TH, A-1101-019, (SKF BEARING: 6001-ZZ DOUBLE SHIELD)	2
B1a — B7	DRVMOTR1205	MOTOR,DRIVE,E610,3700 RPM,RIGHT,MIP,H-1121-002, (MTS-62020078)	1
B1b — B7	DRVMOTR1204	MOTOR,DRIVE,E610,3700 RPM,LEFT,MIP,H-1121-003, (MTS-62010078)	1
B1	INDPART1925	AXLE KEY,0.236"W X 0.236" H X 1.566" L	1
B2	PLSINDV1104	BRAKE,CAP,E610	1
В3	INDPART2056	BRAKE,SOLENOID,E610	1
B4	SCRALNH1037	SCREW,SHCS,M6-1.0 X 16MM	3
B5	SCRALNH1038	SCREW,SHCS,M6-1.0 X 30MM	3
В6	DRVBRSH1006	MOTOR,BRUSH,(MTS: 62400001)	2
В7	PLSINDV1103	MOTOR,BRUSH CAP,(MTS: 62400013)	2
C1a — C11	FRMASMB5479	FREEWHEEL,ASSY,RELEASE LEVER,JAZZY1101/1121,RIGHT	1
C1b — C11	FRMASMB5478	FREEWHEEL,ASSY,RELEASE LEVER,JAZZY1101/1121,LEFT	1
C1a	DWR1101D030	FREEWHEEL,LEVER W/OVERMOLD COVER,RIGHT	1
C1b	DWR1101D031	FREEWHEEL,LEVER W/OVERMOLD COVER,LEFT	1
C2a	DWR1133C014	FREEWHEEL,RELEASE ROD,RIGHT, JAZZY1133,C-1133-014	1
C2b	DWR1133C015	FREEWHEEL,RELEASE ROD,LEFT, JAZZY1133,C-1133-015	1
C3	DWR1114B014	FREEWHEEL,MOTOR RELEASE WELDMENT,JAZZY1113ATS,B-1114-014	1
C4	SCRSETS1023	SCREW,SET,M6-1.0T X 8MM,CUP POINT	1
C5	SPGRAW1046	SPRING,8MM ID X 10MM OD X 10MM TH,COMPRESSION,FREEWHEEL RELEASE LEVER,MTS	1
C6	WASWAVY1008	WASHER,BELLEVILLE,STAINLESS STEEL,0.317"ID X 0.625"OD X .022"TH	1
C7	PINSOTR1000	PIN, COTTER,ZINC PLATED,1/16" X 1"	2
C8	BLTASHL1111	BOLT,SHOULDER,8MM X 45MM,GRADE 12.9,C-1056-025	1
C9	NUTNYLK1019	NUT, NYLOCK,BLACK ZINC,M6-1.0	1
C10	WASNYLN1089	WASHER,NYLON,8MM ID X 19MM OD X 1MM TH	1
C11	WASSTND1100	WASHER,FLAT,ZINC,6MM ID X 14MM OD X 1.5MM TH	1

Notes:

Only for use with the Jazzy1121HD upgrade package. Top-level assemblies will include components designated with corresponding lower-case letters. Parts without lower-case extensions will be included in all top-level assemblies.








REF#	PART#	DESCRIPTION	QTY
A1a — D2	FRMASMB4100	ANTI-TIP,ASSY,CURB CLIMBER FORK,RIGHT,JAZZY1121	1
A1b — D2	FRMASMB4098	ANTI-TIP,ASSY,CURB CLIMBER FORK,LEFT,JAZZY1121	1
A1a — A5	FRMASMB4101	ANTI-TIP,ASSY,BRACKET,CURB CLIMBER, JAZZY1121, RIGHT	1
A1b — A5	FRMASMB4099	ANTI-TIP,ASSY,BRACKET,CURB CLIMBER, JAZZY1121, LEFT	1
A1a	DWR1121D005	CURB CLIMBER WELDMENT,RIGHT, JAZZY1121,D-1121-005	1
A1b	DWR1121D004	CURB CLIMBER WELDMENT,LEFT, JAZZY1121,D-1121-004	1
A2	BLTASHL1087	BOLT,SHOULDER,12MM X 55MM,C-1101-025	1
A3	WASNYLN1080	WASHER,NYLON,13MM ID X 28MM OD X 1.5MM TH	2
A4	SCRHHCS1192	SCREW,HHCS,BLACK ZINC,M8-1.25 X 80MM,GRADE 8.8	2
A5	NUTNYLK1011	NUT,NYLOCK,BLACK ZINC,M8-1.25	2
B1 — B12	FRMASMB4104	SUSPENSION, ASSY, STRUT TOWER, CURB CLIMBER, JAZZY1121	1
B1	DWR1121D007	SUSPENSION,STRUT WELDMENT,1/2-20 X 9.42"L,JAZZY1121,D-1121-007	1
B2	DWR1121B002	SPRING,YELLOW,2.00"L X 1.50"OD X 0.218" WIRE DIA X RATE: 512 LB/,B-1121-002,(LEE SPRING: LHL-1500AB-01)	1
В3	DWR1101B024	PIVOT BUSHING,SUSPENSION STRUT,B-1101-024	1
B4	DWR1121B001	SUSPENSION,SPRING UPPER DISK,JAZZY1121,B-1121-001	1
B5	SCRSHCS1079	SCREW,SHCS,BLACK ZINC,M8-1.25 X 15MM,W/NYLON PATCH	1
В6	DWR1121C002	SPRING,YELLOW,5.00"L X 1.5"OD X 0.206" WIRE DIA X RATE: 95 LBS/S,C-1121-002,(CENTURY SPRINGS: 72863)	1
В7	BLTASHL1085	BOLT,SHOULDER,8MM X 15MM,GRADE 12.9,C-1005-021	2
B8	NUTNJAM1025	NUT,NYLOCK JAM,BLACK ZINC,M10-1.25	1
B9	WASNYLN1085	WASHER,BLACK,NYLON,8.5MM ID X 16MM OD X 3MM TH	2
B10	WASNYLN1084	WASHER,BLACK,NYLON,8.5MM ID X 16MM OD X 1MM TH	2
B11	DWR1121A001	SHIM,BLACK ZINC,1.00" ID X 1.50" OD X 0.25"TH,A-1121-001	1
B12	DWR1121A002	SHIM,BLACK ZINC,0.325" ID X 0.748" OD X 0.14"TH,A-1121-002	1
C1 — C6	FRMASMB4105	WHEEL,ASSY,ANTI-TIP,JAZZY1101/1121	1
C1	DWR1101C013	ANTI-TIP,ADJUSTMENT BRACKET,JAZZY1101/1121,C-1101-013	1
C2	SCRHHCS1222	SCREW,HHCS,BLACK ZINC,M8-1.25 X 75MM,GRADE 10.9	1
C3	NUTNYLK1011	NUT,NYLOCK,BLACK ZINC,M8-1.25	1
C4	SCRHHCS1219	SCREW,HHCS,BLACK ZINC,M6-1.0 X 25MM,GRADE 10.9	1
C5	NUTNYLK1019	NUT,NYLOCK,BLACK ZINC,M6-1.0	2
C6	WASNYLN1085	WASHER,BLACK,NYLON,8.5MM ID X 16MM OD X 3MM TH	1
D1 — D2	WHLASMB1284	WHEEL,ASSY,ANTI-TIP,6",MOLDED,W/BEARING	1
D1	BRGSLV1002	BEARING,BALL,8MM ID X 22MM OD X 7MM TH,608ZZ	2
D2	WHLATIP1001	WHEEL,ANTI-TIP,6",MOLDED	1

Top-level assemblies will include components designated with corresponding lower-case letters. Parts without lower-case extensions will be included in all top-level assemblies.








REF#	PART #	DESCRIPTION	QTY
A1a — D2	FRMASMB6210	ANTI-TIP,ASSY,CURB CLIMBER FORK, JAZZY1121HD, RIGHT	1
A1b — D2	FRMASMB6208	ANTI-TIP,ASSY,CURB CLIMBER FORK, JAZZY1121HD, LEFT	1
A1a — A5	FRMASMB4101	ANTI-TIP,ASSY,BRACKET,CURB CLIMBER, JAZZY1121, RIGHT	1
A1b — A5	FRMASMB4099	ANTI-TIP,ASSY,BRACKET,CURB CLIMBER, JAZZY1121, LEFT	1
A1a	DWR1121D005	CURB CLIMBER WELDMENT,RIGHT, JAZZY1121,D-1121-005	1
A1b	DWR1121D004	CURB CLIMBER WELDMENT,LEFT, JAZZY1121,D-1121-004	1
A2	BLTASHL1087	BOLT,SHOULDER,12MM X 55MM,C-1101-025	1
A3	WASNYLN1080	WASHER,NYLON,13MM ID X 28MM OD X 1.5MM TH	2
A4	SCRHHCS1192	SCREW,HHCS,BLACK ZINC,M8-1.25 X 80MM,GRADE 8.8	2
A5	NUTNYLK1011	NUT,NYLOCK,BLACK ZINC,M8-1.25	2
B1 — B11	FRMASMB6209	ASSY,STRUT,CURB CLIMBER,JAZZY1121HD	1
B1	DWR1121D007	SUSPENSION,STRUT WELDMENT,1/2-20 X 9.42"L,JAZZY1121,D-1121-007	1
B2	DWR1121B002	SPRING,YELLOW,2.00"L X 1.50"OD X 0.218" WIRE DIA X RATE: 512 LB/,B-1121-002,(LEE SPRING: LHL-1500AB-01)	1
В3	DWR1101B024	PIVOT BUSHING, SUSPENSION STRUT, B-1101-024	1
B4	DWR1121B001	SUSPENSION,SPRING UPPER DISK,JAZZY1121,B-1121-001	1
B5	SCRSHCS1079	SCREW,SHCS,BLACK ZINC,M8-1.25 X 15MM,W/NYLON PATCH	1
В6	DWR1121C002	SPRING,YELLOW,5.00"L X 1.5"OD X 0.206" WIRE DIA X RATE: 95 LBS/S,C-1121-002,(CENTURY SPRINGS: 72863)	1
В7	BLTASHL1085	BOLT,SHOULDER,8MM X 15MM,GRADE 12.9,C-1005-021	2
В8	NUTNJAM1025	NUT,NYLOCK JAM,BLACK ZINC,M10-1.25	1
B9	WASNYLN1085	WASHER,BLACK,NYLON,8.5MM ID X 16MM OD X 3MM TH	2
B10	WASNYLN1084	WASHER,BLACK,NYLON,8.5MM ID X 16MM OD X 1MM TH	2
B11	DWR1121A001	SHIM,BLACK ZINC,1.00" ID X 1.50" OD X 0.25"TH,A-1121-001	2
C1 — C6	FRMASMB4105	WHEEL,ASSY,ANTI-TIP,JAZZY1101/1121	1
C1	DWR1101C013	ANTI-TIP,ADJUSTMENT BRACKET,JAZZY1101/1121,C-1101-013	1
C2	SCRHHCS1222	SCREW,HHCS,BLACK ZINC,M8-1.25 X 75MM,GRADE 10.9	1
C3	NUTNYLK1011	NUT,NYLOCK,BLACK ZINC,M8-1.25	1
C4	SCRHHCS1219	SCREW,HHCS,BLACK ZINC,M6-1.0 X 25MM,GRADE 10.9	1
C5	NUTNYLK1019	NUT,NYLOCK,BLACK ZINC,M6-1.0	2
C6	WASNYLN1085	WASHER,BLACK,NYLON,8.5MM ID X 16MM OD X 3MM TH	1
D1 — D2	WHLASMB1284	WHEEL,ASSY,ANTI-TIP,6",MOLDED,W/BEARING	1
D1	BRGSLV1002	BEARING,BALL,8MM ID X 22MM OD X 7MM TH,608ZZ	2
D2	WHLATIP1001	WHEEL,ANTI-TIP,6",MOLDED	1

Top-level assemblies will include components designated with corresponding lower-case letters. Parts without lower-case extensions will be included in all top-level assemblies.





V. Articulating Beam Assembly — Straight




REF#	PART#	DESCRIPTION	QTY
A1 — D3	FRMASMB4107	ARTICULATING BEAM,ASSY,STRAIGHT,W/ CASTER WHEELS, JAZZY1101/1121	1
A1 — A15	FRMASMB4108	ARTICULATING BEAM,ASSY,STRAIGHT,JAZZY1101/1121	1
A1	DWR1101E008	ARTICULATING BEAM,WELDMENT,21"LG,STRAIGHT,JAZZY1101/1121,E-1101-008	1
A2	DWR1101A001	SHIM,BRIGHT ZINC,STEEL,13MM ID X 24MM OD X 2.16MM TH,A-1101-001	2
A3	DWR1101A002	SPACER,BRIGHT ZINC,STEEL,14MM ID X 21MM OD X 23.16MM TH,A-1101-002	2
A4	DWR1101A003	SPACER,BRIGHT ZINC,STEEL,13MM ID X 31.5MM OD X 3.18MM TH,A-1101-003	2
A5	DWR1101A013	O-RING,NEOPRENE,25.4MM ID X 31.8MM OD,A-1101-013	2
A6	DWR1101D008	FRAME,SECOND BEAM CATCH & HANDLE WELDMENT,JAZZY1101/1121,D-1101-008	1
A7	SCRHHCS1220	SCREW,HHCS,BLACK ZINC,M8-1.25 X 45MM,GRADE 10.9	2
A8	DWR1101A012	BEARING,FLANGE,OIL IMPREGNATED BRONZE,12.1MM ID X 19MM OD X 19.1MM TH X 25.4MM FOD X 3.2,A-1101-012	2
A9	DWR1101A011	BEARING,SLEEVE,0.50"ID X 32MM OD X 0.394"TH,A-1101-011,(NGK: 6201-2RS-8)	4
A10	BLTASHL1086	BOLT,SHOULDER,12MM X 100MM	1
A11	NUTNYLK1011	NUT,NYLOCK,BLACK ZINC,M8-1.25	2
A12	DWR1101C014	COVER,ARTICULATING BEAM CAP,ARTICULATING BEAM, JAZZY1101/1121,C-1101-014	2
A13	SCRHHCS1237	SCREW,HHCS,BLACK ZINC,M8-1.25 X 15MM,GRADE 10.9	2
A14	NUTNYLK1073	NUT,NYLOCK,BLACK ZINC,M10-1.5,GRADE 8	1
A15	NUTNYLK1010	NUT,NYLOCK,BLACK ZINC,1/2-20	2
B1a — D3	FRMASMB8621	FORK,ASSY,CASTER,RIGHT,W/GUARD & 8" SOLID CASTER,TEXTURED	1
B1b — D3	FRMASMB8620	FORK,ASSY,CASTER,LEFT,W/GUARD & 8" SOLID CASTER,TEXTURED	1
B1a	DWR1122D010	CASTER COVER,PLASTIC,SMALL,RIGHT,D-1122-010	1
B1b	DWR1122D009	CASTER COVER,PLASTIC,SMALL,LEFT,D-1122-009	1
B2	SCRASBH1014	SCREW,BHSCS,BLACK ZINC,1/4-20 X 0.750"	2
В3	WASSTND1002	WASHER,FLAT,BLACK ZINC,0.3125" ID X 0.6875" OD X 0.0625" TH	2
B4	NUTNJAM1020	NUT,NYLOCK JAM,BLACK ZINC,M8-1.25	1
B5	DWR1010A019	SPACER,STEEL,0.320" ID X 0.875" OD,CUSTOM,A-1010-019	1
B6	WASWAVY1005	WASHER,WAVY,STEEL,0.64" ID X 0.865" OD	1
C1a — C2	FRMASMB8618	FORK,ASSY,CASTER,RIGHT,TEXTURED,C-1122-006	1
C1b — C2	FRMASMB8619	FORK,ASSY,CASTER,LEFT,TEXTURED,C-1122-004	1
C1a	DWR1122C005B	FORK,CASTER,JAZZY1122,RIGHT,TEXTURED,C-1122-005	1
C1b	DWR1122C003B	FORK,CASTER,JAZZY1122,LEFT,TEXTURED,C-1122-003	1
*C2	BRK1206A000	STEM,STEEL,W/1/2-20 NYLOCK JAM NUT,A-1200-006	1
D1 — D3	WHLASMB1388	WHEEL, ASSY, CASTER, 200 X 50 W/6900Z BEARING SLEEVE & SPACER	1
D1	DWR1122A020	SPACER,STEEL,0.432"ID X 0.625"OD X 1.81"TH,13 GAUGE,A-1122-020	1
D2	WHLCSTR1003	WHEEL,CASTER,RIM/TIRE	1
D3	BRGSLV1015	BEARING,BALL,10MM ID X 22MM OD X 6MM TH,6900Z	2

Notes

Top-level assemblies will include components designated with corresponding lower-case letters. Parts without lower-case extensions will be included in all top-level assemblies.

*C2 This item is a quantity of one but is illustrated as two separate parts for assembly purposes.





V. Articulating Beam Assembly — Curved



REF#	PART #	DESCRIPTION	QTY
A1 — D3	FRMASMB6211	ARTICULATING BEAM,ASSY,CURVED,W/ CASTER WHEELS, JAZZY1101/1121	1
A1 — A14	FRMASMB6212	ARTICULATING BEAM,ASSY,CURVED,JAZZY1101/1121	1
A1	DWR1101E030	ARTICULATING BEAM WELDMENT, CURVED, JAZZY1101/1121, E-1101-030	1
A2	DWR1101A001	SHIM,BRIGHT ZINC,STEEL,13MM ID X 25MM OD X 2.5MM TH,A-1101-001	2
A3	DWR1101A002	SPACER,BRIGHT ZINC,STEEL,14MM ID X 21MM OD X 23.16MM TH,A-1101-002	2
A4	DWR1101A003	SPACER,BRIGHT ZINC,STEEL,13MM ID X 31.5MM OD X 3.18MM TH,A-1101-003	2
A5	DWR1101A013	O-RING,NEOPRENE,25.4MM ID X 31.8MM OD,A-1101-013	2
A6	DWR1101D008	FRAME,SECOND BEAM CATCH & HANDLE WELDMENT,JAZZY1101/1121,D-1101-008	1
A7	SCRHHCS1220	SCREW,HHCS,BLACK ZINC,M8-1.25 X 45MM,GRADE 10.9	2
A8	DWR1101A012	BEARING,FLANGE,OIL IMPREGNATED BRONZE,12.1MM ID X 19MM OD X 19.1MM TH X 25.4MM FOD X 3.2,A-1101-012	2
A9	DWR1101A011	BEARING,SLEEVE,0.50"ID X 32MM OD X 0.394"TH,A-1101-011,(NGK: 6201-2RS-8)	4
A10	BLTASHL1086	BOLT,SHOULDER,12MM X 100MM	1
A11	DWR1101C014	COVER,ARTICULATING BEAM CAP,ARTICULATING BEAM, JAZZY1101/1121,C-1101-014	2
A12	SCRHHCS1237	SCREW,HHCS,BLACK ZINC,M8-1.25 X 15MM,GRADE 10.9	2
A13	NUTNYLK1011	NUT, NYLOCK,BLACK ZINC,M8-1.25	2
A14	NUTNYLK1073	NUT,NYLOCK,BLACK ZINC,M10-1.5	1
*B1 — D3	FRMASMB8626	CASTER,ASSY,TEXTURED,W/8"WHEEL,THICK WALL,METRIC,RIGHT,JAZZY1121HD	1
*B1 — D3	FRMASMB8625	CASTER,ASSY,TEXTURED,W/8"WHEEL,THICK WALL,METRIC,LEFT,JAZZY1121HD	1
B1	NUTNJAM1018	NUT,NYLOCK JAM,BLACK ZINC,1/2-20	1
B2	NUTNYLK1011	NUT, NYLOCK,BLACK ZINC,M8-1.25	1
В3	SCRHHCS1261	SCREW,HHCS,BLACK ZINC,M8-1.25 X 90MM,GRADE 10.9	1
B4	WASMETR1010	WASHER,FLAT,BLACK ZINC,M8 X 1.2MM TH	1
C1 — C3	WHLASMB1112	WHEEL,ASSY,CASTER,200 X 50,W608ZZ BEARING SLEEVE/EXTENDED RACE & SPACER	1
C1	IND1222A000	SPACER,BLACK ZINC,STEEL,0.325" ID X 0.63" OD X 1.70" TH,A-1200-022	1
C2	WHLCSTR1003	WHEEL,CASTER,RIM/TIRE	1
C3	BRGSLV1017	BEARING,BALL,8MM ID X 22MM OD X 7MM TH,608ZZ,W/ EXTENDED INNER RACE 9MM	2
D1 — D3	FRMASMB6355	CASTER,ASSY,FORK,TEXTURED,THICK WALL,METRIC W/STEM,NUT,& SPACERS	1
D1	DWR1200C101B	FORK,CASTER,THICK WALL,DOUBLE SIDED,TEXTURED BLACK,C-1200-101	1
*D2	BRK1206A000	STEM,STEEL,W/1/2-20 NYLOCK JAM NUT,A-1200-006	1
D3	SPCROLL1000	SPACER,STEEL,0.3125" ID X 0.4375" OD X 0.125" TH	2


Notes:

This assembly is for use with Jazzy1121HD option and power positioning seating applications.

*B1 — D3 Both sides use identical parts. The two FRMASMB numbers shown above are for installation purposes. Quantities are listed per

*D2 This item is a quantity of one but illustrated as two separate parts for assembly purposes.








REF#	PART#	DESCRIPTION	QTY
A1 — A19	PLSASMB1433	ELECTRONIC,ASSY,VSI TRAY,JAZZY1121,OFFBOARD CHARGER	1
A1	DWR1101E020	ELECTRONIC,TRAY,JAZZY1101/1121	1
A2	DWR1500B045	FUSE MODULE COVER,RUBBER,B-1500-045	1
A3	ELEFUSE1010	FUSE,ATO,7.5AMP	2
A4	DWR1101E007	ELECTRONIC,MOUNTING BRACKET,JAZZY1101/1121	1
A5	SCRPPH1316	SCREW,PPH,BLACK ZINC,M58 X 15MM,(MTS: 50105036)	2
A6	ELEASMB2869	ELECTRONIC,ASSY,BATTERY CABLE,BLACK,W/125A/32V FUSEABLE LINK,QUICK DISCONNECT.H-1122-022	2
A7	DWR1101C038	AC CAP,BELLY PAN SIDE,SNAP-IN,C-1101-038	2
A8	SCRPPH1266	SCREW,PPH,ZINC,M35 X 16MM	4
A9	DWR1000B297	CAP,PLASTIC,BLANK,3.32" X 1.60",JAZZY ELECTRONICS TRAY,B-1000-297	1
A10	HARUSHD2080	ELECTRONIC,HARNESS,3 PIN,RUN PLUG,FEMALE,H-1200-071	1
A11	DWR1010H021	ELECTRONIC,HARNESS,POWER INTERFACE,PILOT/VSI,H-1101-021	1
A12	SCRPPH1044	SCREW,PPH,BLACK ZINC,M47 X 12MM	6
A13	PLSTIE1004	TIE WRAP,BLACK,8"	9
A14	SCRPPH1148	SCREW,PPH,BLACK ZINC,M35 X 10MM	6
A15	ELEASMB1000	ELECTRONIC,ASSY,CIRCUIT BREAKER,50 AMP,COMPLETE	1
A16	NUTNYLK1072	NUT,NYLOCK,BLACK ZINC,M35	6
A17	WASMETR1094	WASHER,FLAT,BLACK ZINC,5MM ID X 11MM OD X 1.5MM TH	2
A18	ELEASMB2892	ELECTRONIC,ASSY,FUSE MODULE,W/CONTACTS,QUANTUM BLAST,H-1500-016	1
A19	DWR1000C245	MOUNTING BRACKET,BEAU CONNECTOR MOUNTING BRACKET,C-1000-245	1





VI. Utility Tray Assembly — Remote Plus




REF#	PART #	DESCRIPTION	QTY
A1 — B19	ELEASM B3153	ELECTRONIC, ASSY, REM OTE PLUS, JAZZY1121, OFFBOARD CHARGER	1
A1	DWR1010H025	ELECTRONIC,HARNESS,BATTERY, REMOTE PLUS, JAZZY1101/1121	1
A2	DWR1010H026	ELECTRONIC, HARNESS, MOTOR CABLE, REMOTE PLUS, RIGHT, JAZZY1101/1121	1
A3	DWR1010H027	ELECTRONIC, HARNESS, MOTOR CABLE, REMOTE PLUS, LEFT, JAZZY1101/1121	1
A4	ELEASMB2869	ELECTRONIC, ASSY, BATTERY CABLE, W/125A/32V FUSEABLE LINK, QUICK DISCONNECT, H-1122-022	2
A5	PLSTIE1004	TIE WRAP,BLACK,8"	2
A6	SCRASBH1080	SCREW,BHSCS,BLACK ZINC,M47 X 30MM	2
B1 — B19	PLSASMB1448	ELECTRONIC, ASSY, REMOTE PLUS TRAY, OFFBOARD CHARGER, JAZZY1121INT	1
B1	DWR1101E020	ELECTRONIC,TRAY,JAZZY1101/1121	1
B2	ELEASMB1000	ELECTRONIC, ASSY, CIRCUIT BREAKER, 50 AMP, COMPLETE	1
В3	DWR1500B045	FUSE MODULE COVER,RUBBER,B-1500-045	1
B4	ELEASMB2892	ELECTRONIC, ASSY, FUSE MODULE, W/CONTACTS, QUANTUM BLAST, H-1500-016	1
B5	ELEFUSE1010	FUSE,ATO,7.5AMP	2
B6	DWR1000B297	CAP,PLASTIC,BLANK,3.32" X 1.60",JAZZY ELECTRONICS TRAY,B-1000-297	1
B7	DWR1101E007	ELECTRONIC,MOUNTING BRACKET,JAZZY1101/1121	1
B8	DWR1101C016	CAP,PLASTIC,BEAU CONNECTOR/REMOTE CAP,C-1101-016	1
В9	SCRPPH1148	SCREW,PPH,BLACK ZINC,M35 X 10MM	6
B10	NUTNYLK1072	NUT,NYLOCK,BLACK ZINC,M35	6
B11	DWR1101C038	AC CAP,BELLY PAN SIDE,SNAP-IN,C-1101-038	2
B12	DWR1010H051	ELECTRONIC, HARNESS, POWER INTERFACE/CHARGER CABLE, REMOTE PLUS, RIGHT, H-1010-051	1
B13	PLSTIE1004	TIE WRAP,BLACK,8"	9
B14	SCRPPH1316	SCREW,PPH,BLACK ZINC,M58 X 15MM,(MTS: 50105036)	2
B15	SCRPPH1044	SCREW,PPH,BLACK ZINC,M47 X 12MM	4
B16	SCRPPH1032	SCREW,PPH-THREAD CUTTING,BLACK ZINC,6-32 X 0.375",TYPE F	4
B17	SCRPPH1266	SCREW,PPH,ZINC,M35 X 16MM	4
B18	WASMETR1094	WASHER,FLAT,BLACK ZINC,5MM ID X 11MM OD X 1.5MM TH	2
B19	HARUSHD2080	ELECTRONIC,HARNESS,3 PIN,RUN PLUG,FEMALE,H-1200-071	1
*C1a	ELEASM B5044	CONTROLLER, REMOTE PLUS, ASSY, 80 AMP, JAZZY SERIES PROGRAM	1
*C1b	ELEASM B5045	CONTROLLER, REMOTE PLUS (PROFILES), ASSY, 80 AMP, JAZZY SERIES PROGRAM	1

Notes:

*C1a Used with single profile controllers.*C1b Used with multiple profile controllers.








REF#	PART#	DESCRIPTION	QTY
A1	CTLDC1270	CONTROLLER,VSI,4 KEY,50 (70)AMP,W/RIGHT ANGLE CONNECTOR,H-1010-040,(P&G: D50148)	1
A2	ELEASMB4956	CONTROLLER,VSI,ASSY,4 KEY,70 AMP,W/RIGHT ANGLE CONNECTOR,JAZZY1121HD PROGRAM	1
B1	CTLDC1271	CONTROLLER,VSI,6 KEY,50 (70)AMP,W/RIGHT ANGLE CONNECTOR,H-1010-041,(P&G: D50149)	1
B2	CTLDC1382	CONTROLLER,VSI,6 KEY,70 AMP,1 ACTUATOR,W/RIGHT ANGLE CONNECTOR,H-1010-081,(P&G: D50530)	1
C1 — D2	FRMASMB4767	JOYSTICK MOUNTING BRACKET, ASSY, VSI CONTROLLER, INLINE-RIGHT, C-1010-193	1
C1	SCRPPH1122	SCREW,PPH,BLACK ZINC,M58 X 18MM	2
D1 — D2	FRMASMB6220	JOYSTICK MOUNTING BRACKET, ASSY, VSI CONTROLLER, INLINE-RIGHT, W/END CAPS, E-1010-280	1
D1	DWR1010D123	MOUNTING BRACKET, JOYSTICK WELDMENT, VSI, INLINE RIGHT, D-1010-123	1
D2	PROECAP1014	CAP, END,PLASTIC,3/4" X 16GA,ROUND	2
E1	HARUSHD2055	ELECTRONIC,HARNESS,1M LONG,DX BUS,DYNAMIC	1
E2	CTLDC1032	CONTROLLER, DYNAMIC, JSTK, DX-REMG805	1
F1 — F2	FRMASMB1789	FIXED INLINE JOYSTICK MOUNT, ASSY, EUROPA	1
F1	DWG1200B270	MOUNTING BRACKET, CONTROLLER, DYNAMIC DX JOYSTICK, B-1200-027	1
F2	SCRPFH1005	SCREW, PFH,BLACK ZINC,10-32 X .0375"	4
G1	CTLDC1112	CONTROLLER, JOYSTICK MODULE, REMOTE PLUS, 3 KEY, SINGLE PROFILE, (P & G: D49671)	1
G2	CTLDC1368	CONTROLLER, JOYSTICK MODULE, REMOTE PLUS, 3 KEY W/PROFILES, FLYING LEAD, (P&G D50524)	1
H1 — H2	FRMASMB2066	FIXED INLINE JOYSTICK MOUNT, ASSY, REMOTE PLUS	1
H1	DWR1200C103	MOUNTING BRACKET, JOYSTICK WELDMENT, REMOTE PLUS, INLINE, C-1200-103	1
H2	SCRPPH1006	SCREW,PPH,8-32 X 1/2",G2,ELMR	2
I1	HARUSHD2189	ELECTRONIC,HARNESS,1.16 LONG,REDEL CABLE, REMOTE PLUS, JOYSTICK TO TRAY	1
J1	CTLDC1369	CONTROLLER,MICRODRIVE REMOTE PLUS JOYSTICK MODULE,(P&G: D50470)	1
K1	PLSINDV1164	VSI SHIELD,(P&G: P76630)	1

Refer to the Quantum Accessories IPB for additional Remote Plus and Europa Controllers.








REF#	PART#	DESCRIPTION	QTY
A1 — B5	FRMASMB4112	FOOT PLATFORM,ASSY,GRID STYLE,9.55"D,W/MOUNTING BRACKET,JAZZY1101/1121	1
A1 — A12	FRMASMB4113	MOUNTING BRACKET, ASSY, FOOT PLATFORM, JAZZY1101/1121	1
A1	DWR1101C010	FOOT PLATFORM MOUNTING BRACKET, JAZZY1101/1121	1
A2	SCRHHCS1258	SCREW,HHCS,BLACK ZINC,M6-1.0 X 40MM,GRADE 10.9	2
A3	NUTNJAM1021	NUT,NYLOCK JAM,BLACK ZINC,M6-1.0	2
A4	PINCLEV1031	QUICK RELEASE PIN,6MM DIA X 38.1MM GRIP LENGTH,W/"T" HANDLE,B-1101-014	1
A5	DWR1010B090	SPACER, STEEL,0.328" ID X 0.50" OD X 1.00" TH,B-1010-090	1
A6	WASNYLN1066	WASHER,BLACK,NYLON,8MM ID X 22MM OD X 3MM TH	1
A7	DWR1000B458	FOOT PLATFORM MOUNTING SLEEVE,FLANGE,STEEL,B-1000-458	1
A8	DWR1000B459	SPACER,STEEL,0.328" ID X 0.50" OD X 0.076" TH,B-1000-459	1
A9	SCRSHCS1170	SCREW,SHCS,BLACK ZINC,M8-1.25 X 80MM,GRADE 12.9	1
A10	PROECAP1080	CAP,VINYL,0.50" ID,FITS 0.50" - 0.5625" OD,0.50" H,(MCMASTER CARR: 9753K44)	2
A11	NUTNJAM1020	NUT,NYLOCK JAM,BLACK ZINC,M8-1.25	1
*A12	BRGSBTH1000	BEARING,BALL THRUST,STAINLESS STEEL,0.3125" ID X 0.625" OD X 0.193" TH,(MCMASTER-CARR: 6655K34)	1
B1 — B5	FRMASMB2295	FOOT PLATFORM,ASSY,GRID STYLE,9.55"D,W/MAT & HARDWARE,E-1000-108	1
B1	DWR1000E106	FOOT PLATFORM BASE,ALUMINUM,9.55"D,GRID STYLE,E-1000-106	1
B2	DWR1000E107	MAT,FOOT PLATFORM,RUBBER,9.55"D,GRID STYLE,E-1000-107	1
В3	SCRPPH1080	SCREW,PPH,BLACK ZINC,8-32 X 0.75"	7
B4	SCRSETS1090	SCREW, SET,BLACK ZINC,3/8-24 X 1.50",W/NYLON PATCH	1
B5	WASSAE1008	WASHER,SAE,BLACK ZINC,#8,0.1875" ID X 0.4375" OD X 0.0468" TH	7

*A12 This item is a quantity of one but is illustrated as three parts for assembly purposes.








REF#	PART#	DESCRIPTION	QTY
A1 — B5	FRMASMB4111	FOOT PLATFORM,ASSY,GRID STYLE,11.62"D,W/MOUNTING BRACKET,JAZZY1101/1121	1
A1 — A12	FRMASMB4113	MOUNTING BRACKET, ASSY, FOOT PLATFORM, JAZZY1101/1121	1
A1	DWR1101C010	FOOT PLATFORM MOUNTING BRACKET,C-1101-010	1
A2	SCRHHCS1258	SCREW,HHCS,BLACK ZINC,M6-1.0 X 40MM,GRADE 10.9	2
A3	NUTNJAM1021	NUT,NYLOCK JAM,BLACK ZINC,M6-1.0	2
A4	PINCLEV1031	QUICK RELEASE PIN,6MM DIA X 38.1MM GRIP LENGTH,W/"T" HANDLE,B-1101-014	1
A5	DWR1010B090	SPACER, STEEL,0.328" ID X 0.50" OD X 1.00" TH,B-1010-090	1
A6	WASNYLN1066	WASHER,BLACK,NYLON,8MM ID X 22MM OD X 3MM TH	1
A7	DWR1000B458	FOOT PLATFORM MOUNTING SLEEVE,FLANGE,STEEL,B-1000-458	1
A8	DWR1000B459	SPACER,STEEL,0.328" ID X 0.50" OD X 0.076" TH,B-1000-459	1
A9	SCRSHCS1170	SCREW,SHCS,BLACK ZINC,M8-1.25 X 80MM,GRADE 12.9	1
A10	PROECAP1080	CAP,VINYL,0.50" ID,FITS 0.50" - 0.5625" OD,0.50" H,(MCMASTER CARR: 9753K44)	2
A11	NUTNJAM1020	NUT,NYLOCK JAM,BLACK ZINC,M8-1.25	1
*A12	BRGSBTH1000	BEARING,BALL THRUST,STAINLESS STEEL,0.3125" ID X 0.625" OD X 0.193" TH,(MCMASTER-CARR: 6655K34)	1
B1 — B5	FRMASMB2117	FOOT PLATFORM, ASSY, GRID STYLE, 11.62"D, W/MAT & HARDWARE, E-1000-105	1
B1	DWR1000E103	MAT,FOOT PLATFORM,11.62"D,GRID STYLE,E-1000-103	1
B2	DWR1000E104	FOOT PLATFORM BASE,ALUMINUM,11.62"D,GRID STYLE,E-1000-104	1
В3	SCRPPH1080	SCREW,PPH,BLACK ZINC,8-32 X 0.75"	8
B4	SCRSETS1090	SCREW, SET,BLACK ZINC,3/8-24 X 1.50",W/NYLON PATCH	1
B5	WASSAE1008	WASHER,SAE,BLACK ZINC,#8,0.1875" ID X 0.4375" OD X 0.0468" TH	8

*A12 This item is a quantity of one but is illustrated as three parts for assembly purposes.








REF#	PART#	DESCRIPTION	QTY
A1 — B3	FRMASMB7888	FOOT PLATFORM,ASSY,HD STEEL,13.5"W,11"D,W/MOUNTING BRACKET,JAZZY1121HD	1
A1 — A12	FRMASMB4113	MOUNTING BRACKET,ASSY,FOOT PLATFORM,JAZZY1101/1121	1
A1	DWR1101C010	FOOT PLATFORM MOUNTING BRACKET,C-1101-010	1
A2	NUTNJAM1021	NUT,NYLOCK JAM,BLACK ZINC,M6-1.0	2
A3	PINCLEV1031	QUICK RELEASE PIN,6MM DIA X 38.1MM GRIP LENGTH,W/"T" HANDLE,B-1101-014	1
A4	SCRHHCS1258	SCREW,HHCS,BLACK ZINC,M6-1.0 X 40MM,GRADE 10.9	2
A5	WASNYLN1066	WASHER,BLACK,NYLON,8MM ID X 22MM OD X 3MM TH	1
A6	NUTNJAM1020	NUT,NYLOCK JAM,BLACK ZINC,M8-1.25	1
*A7	BRGSBTH1000	BEARING,BALL THRUST,STAINLESS STEEL,0.3125" ID X 0.625" OD X 0.193" TH,(MCMASTER-CARR: 6655K34)	1
A8	DWR1010B090	SPACER, STEEL,0.328" ID X 0.50" OD X 1.00" TH,B-1010-090	1
A9	DWR1000B458	FOOT PLATFORM MOUNTING SLEEVE,FLANGE,STEEL,B-1000-458	1
A10	DWR1000B459	SPACER,STEEL,0.328" ID X 0.50" OD X 0.076" TH,B-1000-459	1
A11	SCRSHCS1170	SCREW,SHCS,BLACK ZINC,M8-1.25 X 80MM,GRADE 12.9	1
A12	PROECAP1080	CAP,VINYL,0.50" ID,FITS 0.50" - 0.5625" OD,0.50" H,(MCMASTER CARR: 9753K44)	2
B1 — B3	FRMASMB7887	FOOT PLATFORM,ASSY,HD STEEL,13"W,11"D,W/MAT & HARDWARE,E-1000-191	1
B1	DWR1000E183	FOOT PLATFORM WELDMENT,HD,LARGE,E-1000-183	1
B2	DWR1000C376	MAT,RUBBER,BLACK,VINYL,V-RIB,C-1000-376	1
В3	SCRSETS1090	SCREW, SET,BLACK ZINC,3/8-24 X 1.50",W/NYLON PATCH	1


*A7 This item is a quantity of one but is illustrated as three parts for assembly purposes.






REF#	PART #	DESCRIPTION	QTY
A1 — B5	WHLASMB1475	WHEEL,ASSY,FLAT FREE DRIVE WHEEL,3.00 X 8",HONEYCOMB,3 PIECE RIM W/FORGED STEEL 17MM KEYWAY,C-9210-1	2
A1	TIRSOLI1036	TIRE,FLAT FREE,300-8 X 14",DRIVE WHEEL, AQUA-TREAD,C9210	1
B1 — B5	WHLASMB1502	WHEEL,ASSY,RIM,HONEYCOMB, ALUMINUM,8",HONEYCOMB	1
B1	DWR1200E221	WHEEL,RIM,ALUMINUM,3.00-8",OUTER, 8" HONEYCOMB, W/HUB MOUNTING HOLES,E-1200-221	1
B2	DWR1200D234	WHEEL,HUB,FORGED STEEL,2.756"BC X 3.386"DIA X 1.89"TH,17MM BORE & KEYWAY, FORGED STEEL, HONEYCOMB,D-1200-234	1
В3	WHLRIMS1051	WHEEL,RIM,ALUMINUM MAG,JAZZY1120/1170	1
B4	SCRSHCS1088	SCREW,SHCS,BLACK OXIDE,5/16-24 X1.00",W/OIL RESISTANT LOCTITE	4
B5	SCRALNH1033	SCREW,SHCS,BLACK ZINC,1/4-20 X 0.75",W/NYLON PATCH	8





REF#	PART#	DESCRIPTION	QTY
A1 — B5	WHLASMB1476	WHEEL, ASSY, PNEUMATIC DRIVE WHEEL, 3.00 X 8", HONEYCOMB, 3-PIECE RIM W/FORGED STEEL 17MM KEYWAY, C-9201-1	2
A1	TUBSTND1002	TUBE,300-8,ANGLED VALVE W/CAP	1
A2	TIRPNEU1032	TIRE,PNEUMATIC,300-8 X 14",DRIVE WHEEL, AQUA-TREAD,C9210	1
A3	DWR7700A010	SLIP COVER,VALVE STEM,35PSI,2.4 BAR,MAX RATING,A-7700-010	1
B1 — B5	WHLASMB1502	WHEEL,ASSY,RIM,HONEYCOMB, ALUMINUM,8",HONEYCOMB	1
B1	WHLRIMS1051	WHEEL,RIM,ALUMINUM MAG,JAZZY1120/1170	1
B2	DWR1200E221	WHEEL,RIM,ALUMINUM,3.00-8",OUTER, 8" HONEYCOMB, W/HUB MOUNTING HOLES,E-1200-221	1
В3	DWR1200D234	WHEEL,HUB,FORGED STEEL,2.756"BC X 3.386"DIA X 1.89"TH,17MM BORE & KEYWAY, FORGED STEEL, HONEYCOMB,D-1200-234	1
B4	SCRSHCS1088	SCREW,SHCS,BLACK OXIDE,5/16-24 X1.00",W/OIL RESISTANT LOCTITE	4
B5	SCRALNH1033	SCREW,SHCS,BLACK ZINC,1/4-20 X 0.75",W/NYLON PATCH	8






REF#	PART #	DESCRIPTION	QTY
A1a — A4	SHDASM B3429	SHROUD,ASSY,RED,JAZZY1101/1121	1
A1b — A4	SHDASM B3430	SHROUD,ASSY,BLUE,JAZZY1101/1121	1
A1c — A4	SHDASM B3435	SHROUD,ASSY,CUSTOM COLOR,JAZZY1101/1121	1
A1a	SHDPNTD2056	SHROUD,PAINTED,RED,JAZZY1101/1121	1
A1b	SHDPNTD2057	SHROUD,PAINTED,BLUE,JAZZY1101/1121	1
A1c	SHDPNTD2062	SHROUD,PAINTED,CUSTOM COLOR,JAZZY1101/1121	1
A2	SCRPPSM1011	SCREW,PPH-SHEET METAL,BLACK ZINC,#8 X 0.50"	4
А3	DWR1101E018	SHROUD,CENTER,JAZZY1101/1121	1
A4	CLISTL1013	CLIP,STUD,1"L,STRAIGHT	4

Custom color available. Please contact Pride Technical Service for more information. Top-level shroud assemblies contain illustrated painted shrouds designated with corresponding low er-case extensions. Painted shrouds are not available individually. The assembly must be purchased. Parts without low er-case extensions will be included in all top-level assemblies.






REF#	PART #	DESCRIPTION	QTY
A1a — A5	SHDASM B3439	SHROUD,ASSY,RED,JAZZY1101/1121, POWER SEAT OPTION	1
A1b — A5	SHDASMB3440	SHROUD,ASSY,BLUE,JAZZY1101/1121, POWER SEAT OPTION	1
A1c — A5	SHDASMB3445	SHROUD, ASSY, CUSTOM COLOR, JAZZY1101/1121, POWER SEAT OPTION	1
A1a	SHDPNTD2056	SHROUD,PAINTED,RED,JAZZY1101/1121	1
A1b	SHDPNTD2057	SHROUD,PAINTED,BLUE,JAZZY1101/1121	1
A1c	SHDPNTD2062	SHROUD,PAINTED,CUSTOM COLOR,JAZZY1101/1121	1
A2	SCRPPSM1011	SCREW,PPH-SHEET METAL,BLACK ZINC,#8 X 0.50"	2
А3	DWR1101D021	SHROUD,CENTER,W/POWER SEAT CUT-OUT, JAZZY1101/1121	1
A4	DWR9500C013	MAT,RUBBER,SCUFF GUARD,BLACK,SEAT POST, FLANGE,C-9500-013	1
A5	CLISTL1013	CLIP,STUD,1"L,STRAIGHT	4

Custom color option available. Please contact Pride Technical Service for more information. Top-level shroud assemblies contain illustrated painted shrouds designated with corresponding low er-case extensions. Painted shrouds are not available individually. The assembly must be purchased. Parts without low er-case extensions will be included in all top-level assemblies.





REF#	PART#	DESCRIPTION	QTY
A1	ELECHG1010	CHARGER,OFFBOARD,110V-230V,5 AMP,PRESET TO 230V,H-1010-114,(MTS:HP-8204B)	1
B1	ELEFUSE1020	FUSE,FAST ACTING,10AMP,250V,5MM OD X 20MM L,(MULTICOMP: 92N5130)	1

XII. Owner's and Decal Packages



Owner's Package

REF#	PART #	DESCRIPTION	QTY
*N/A	INDASM B1256	OWNER'S PACKAGE,ASSY,JAZZY1121INT, UNITED KINGDOM	1
*N/A	PROZBAG1003	BAG,10"L X 12"W X 2MM,ZIPLOCK	1
*N/A	HARUSHD2041	ELECTRONIC,HARNESS, POWER CORD,AC W/UK PLUG,H-1020-049	1
*N/A	INFM ANU1935	MANUAL,OWNER'S,JAZZY1121, UNITED KINGDOM	1
*N/A	VELDBLS1001	REUSABLE FASTENER,1/2"W X 8"L,HOOK & LOOP	1
*N/A	INFINFB1246	CONTROLLER, CONTROLLER ROUTING INSTALLATION INSTRUCTIONS, JAZZY1101/1121	1
*N/A	DWR1020L014	HANG TAG,ROUND,"MADE IN TAIWAN", DISTRIBUTED GLOBALLY,L-1020-014	1
*N/A	PLSTIE1013	CABLE TIE,BLACK,NYLON,8" L,40 LBS TENSILE STRENGTH	5
*N/A	INFINFB1148	MANUAL, WARRANTY CARD, UNITED KINGDOM	1
*N/A	DWR1020L003	HANG TAG, FULLY RECHARGE BATTERY BEFORE INITIAL USE, 6-1/4" X 3-3/8", L-1020-003	1
*N/A	INFINFB1841	MANUAL, INSTALLATION INSTRUCTIONS, CONTROLLER CABLE ROUTING, CONTOURED SEAT	1


Standard Decal Package

REF#	PART#	DESCRIPTION	QTY
*N/A	DECASMB1128	DECAL,ASSY,JAZZY1121INT, UNITED KINGDOM	1
*N/A	DECINDV1041	DECAL,CE,ROUND,5/8",WHITE/BLACK,L-1234-190	1
*N/A	DWG1000L010	DECAL, CAUTION - OPERATIONAL INSTRUCTION, ARM COVER	1
*N/A	DWR1122L004	DECAL,BATTERY CONNECTION, JAZZY1122	1
*N/A	DWR1000L054	DECAL,WARNING,MAXIMUM WEIGHT CAPACITY, 300LBS/135KG	1
*N/A	DWR1000L033	DECAL,"! READ OWNER'S MANUAL BEFORE OPERATING POWERCHAIR" ENGLISH	1
*N/A	DWR1050L005	DECAL, POWERCHAIR, US AND FOREIGN PATENTS PENDING, MFG IN TAIWAN	1
*N/A	DWR1000L034	DECAL,FREEWHEEL WARNING,ENGLISH	1
*N/A	DWR1121L001	DECAL,NAME, JAZZY1121, RIGHT SIDE	1
*N/A	DWR1121L002	DECAL,NAME, JAZZY1121, LEFT SIDE	1
*N/A	DWR1121L003	DECAL,PRIDE MOBILITY PRODUCTS,WEBSITE,L-1121-003	1
*N/A	DWR1000L129	DECAL, WARNING, SEALED "LEAD ACID (SLA) OR GEL-CEL BATTERIES ONLY"	1
*N/A	DWR1234L624	DECAL,ICON,WARNING,CONTROLLER CORD PLACEMENT,ENGLISH,L-1234-624	1
*N/A	DWR1234L957	DECAL,WARNING,DEFLATE TIRE COMPLETELY,L-1234-957	2


HD Decal Package

REF#	PART #	DESCRIPTION	QTY
A1- A10	DECASMB1340	DECAL PACKAGE, JAZZY1121HD, DOMESTIC	1
N/A	DWG1000L010	DECAL CAUTION - OPERATIONAL INSTRUCTION, ARM COVER	1
N/A	DWR1122L004	DECAL, BATTERY CONNECTION, JAZZY1122	1
N/A	DWR1000L057	DECAL, WARNING, MAXIMUM WEIGHT CAPACITY, 400 LBS./180 KG	1
N/A	DWR1000L033	DECAL, "! READ OWNER'S MANUAL BEFORE OPERATING POWERCHAIR" ENGLISH	1
N/A	DWR1050L005	DECAL, POWERCHAIR, US AND FOREIGN PATENTS PENDING, MFG IN TAIWAN	1
N/A	DWR1000L034	DECAL,FREEWHEEL WARNING,ENGLISH	1
N/A	DWR1000L014	DECAL, EMI-RFI, 20V/M, 2001	1
N/A	DWR1121L001	DECAL, NAME, JAZZY1121, RIGHT SIDE	1
N/A	DWR1121L002	DECAL, NAME, JAZZY1121, LEFT SIDE	1
N/A	DWR1121L003	DECAL, PRIDE LOGO W/ WEBSITE, L-1121-003	1











Pride Mobility Products® Corporation

182 Susquehanna Avenue Exeter, PA 18643 USA

Pride Mobility Products® Company

380 Vansickle Road, Unit 350 St. Catharines, Ontario L2R 6P7 Canada

Pride Mobility Products® Ltd.

32 Wedgwood Road Bicester, Oxon OX26 4UL UK

Pride Mobility Products® Australia Pty. Ltd.

21 Healey Road Dandenong, 3175 Victoria, Australia

Pride Mobility Products® Italia S.r.I.

Via del Porgresso - ang. Via del Lavoro Loc. Prato della Corte 00065 - Fiano Romano (RM)

Pride Mobility Products® Europe B.V.

Castricummer Werf 26 1901 RW Castricum The Netherlands